

ANASTOMOSIS LEAK TESTER-OP

Black for OPEN surgery
inflates to ≈ 30 cmH₂O

Catalog No. ALT-30



ANASTOMOSIS LEAK TESTER-LP

White for LAPAROSCOPIC surgery
inflates to ≈ 30 cmH₂O above
intraoperative pressure when
abdomen is insufflated to 15 mmHg

Catalog No. ALT-50



U.S. Patent No. 8,216,159

- Simplify rectal inflation
- Completely occlude the anus
- Signal when targeted pressure is reached
- Prevent overdistention
- Compensate for loss of tactile sensation in laparoscopic surgery
- Sterile, compact and disposable



8 Technology Drive
East Setauket, NY 11733
P: 631-444-4448
F: 631-444-5330
orders@toolsforsurgery.com

www.toolsforsurgery.com

Colorectal Anastomosis Leak Testers

**EXPEDITE AND
STANDARDIZE
THE 'BUBBLE TEST'**

**CONTROL AND DOCUMENT
INFLATION PRESSURE**



ANASTOMOSIS LEAK TESTER-LP

Air leak testing is recommended for all left-sided colorectal anastomoses.*



www.toolsforsurgery.com/ALT

* Ricciardi R et al. Anastomotic leak testing after colorectal resection: what are the data? Arch Surg 2009;144(5):407-411.
Kwon S et al. Routine leak testing in colorectal surgery in the Surgical Care and Outcomes Assessment Program. Arch Surg 2012;147(4):345-351.